easy-way-coil[®]
spring system

Leading in Orthodontic Implants

psm-medical.com

Table of Content | Inhaltsverzeichnis

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


Predictable | Berechenbar

Introduction | Einführung

The **easy-way-coil® spring system** is a highly effective device with a stainless steel spring at its centre. This ensures an extremely long activation distance, offers different directions of force and easy re-activation. A clinically proven technique for the effective alignment of displaced incisors, canines and premolars.

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
easy-way-coil® spring system | easy-way-coil® Spring System

Multi-Directional – Precise Tooth Movement in any desired Direction

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Activation Distance = 2 mm Optimum Force = 32 cN

- Precisely quantified force amounts
- Intermittent force delivery
- Variable direction of movement
- Unilateral or bilateral
- Use with lingual and buccal brackets
- Simple reactivation
- Clinically reliable & economical
- High level of comfort for the patient

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Passive retaining spring | Passive Platzhalterfeder

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Before using the tension spring, you must ensure that there is sufficient space for the tooth to be inserted. With the .018" method for the bracket slot, a stainless steel arch wire of .016 x .022" is sufficient for anchoring, corresponding to a .020 x .020" arch wire of the same quality for the .022" method. In the area of the prepared gap, a precisely fitting piece of the spring is pushed onto the anchoring arch as a placeholder. The first activation takes place seven days after the operation in 4 steps (see Fig. 1-4). Once the canine has broken through and the spring has been shortened to such an extent that it can no longer be used (see Fig. 5), it is carefully removed from the lingual button after loosening the ligature wire by gently pulling it laterally with How-pliers or similar. If the button breaks off in the process, a new button can be glued in the same place. The further procedure in Fig. 6 is shown by a corresponding C-chain, which leads from the bracket of the 2nd premolar palatally to the buttons of the canine and on to the bracket of the 1st premolar. This process is used for further buccal movement of the canine.

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The patented stainless steel spring of the **EWC Spring system** guarantees a precisely adjustable force level with discontinuous force duration. The tension spring is shortened every 4 weeks for reactivation. 1 mm shortening (corresponds to 3 spiral rings on the spring) generates a force of approx. 16 cN, 2 mm shortening corresponds to approx. 32 cN, depending on the required force.

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PSM Medical GmbH
PSM North America, Inc.